

Product Datasheet

Cytochrome C Blocking Peptide (orb706529)

Description The peptide is used to block Anti-Cytochrome c antibody

Reactivity Human, Mouse, Porcine, Rat, Sheep

Tested Applications Blocking

Target CYCS

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID **P99999, P62897, P62898**

Entrez **54205, 100363502, 13063, 25309**

Expiration Date 12 months from date of receipt.