

Product Datasheet

ASK1 (phospho-T838) Blocking Peptide (orb706553)

Description	The peptide is used to block Anti-ASK1 (pT838) antibody <a href= https://www.biorbyt.com/ask1-pt838-a...
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	MAP3K5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99683, O35099
Entrez	4217, 26408
Expiration Date	12 months from date of receipt.