

Product Datasheet

FHIT (phospho-Y114) Blocking Peptide (orb706572)

Description	The peptide is used to block Anti-FHIT (pY114) antibody
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	FHIT
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O89106, P49789
Entrez	14198, 2272
Expiration Date	12 months from date of receipt.