

Product Datasheet

GABBR2 (phospho-S893) Blocking Peptide (orb706576)

Description The peptide is used to block Anti-GABBR2 (pS893) antibody

Reactivity Human, Mouse, Rat

Tested Applications Blocking

Target GABBR2

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID **Q80T41, O88871, O75899**

Entrez **242425, 83633, 9568**

Expiration Date 12 months from date of receipt.