

Product Datasheet

SIX5 Blocking Peptide (orb706675)

Description The peptide is used to block Anti-SIX5 antibody
<a href=https://www.biorbyt.com/six5-antibody-orb665...

Reactivity Human, Mouse

Tested Applications Blocking

Target SIX5

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID **P70178, Q8N196**

Entrez **147912, 20475**

Expiration Date 12 months from date of receipt.