

Product Datasheet

Frizzled 5 Blocking Peptide (orb706705)

Description The peptide is used to block Anti-Frizzled 5 antibody

Reactivity Human, Mouse, Rat

Tested Applications Blocking

Target FZD5; FZD8

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID [Q8CHL0](#), [Q13467](#), [Q498S8](#), [Q9H461](#), [Q9EQD0](#), [Q61091](#)

Entrez [14370](#), [8325](#), [7855](#), [100909849](#), [14367](#), [317674](#)

Expiration Date 12 months from date of receipt.