

Product Datasheet

IgA1/2 Blocking Peptide (orb706721)

Description The peptide is used to block Anti-IgA1 antibody
<a href=https://www.biorbyt.com/iga1-2-antibody-orb6...

Reactivity Human

Tested Applications Blocking

Target IGHA1; IGHA2

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID **P01877, P01876**

Entrez **3493, 3494**

Expiration Date 12 months from date of receipt.