

Product Datasheet

Actin-pan (AcK52/50/51) Blocking Peptide (orb706734)

Description	The peptide is used to block Anti-Actin-pan (AcK52/50/51) antibody <a href=https://www.biorbyt.com/a...
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	ACTA2; ACTB; ACTC1; ACTG1; ACTG2; ACTA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62736, P62738, P60711, P63259, P68032, P68033, P68035, P63261, P63267, P62737, P63260, P60710, P68133, P60709
Entrez	72, 11475, 71, 59, 11464, 29275, 100361457, 60, 11465, 81633, 11461, 70, 81822, 58
Expiration Date	12 months from date of receipt.