

Product Datasheet

Human THBD protein (orb706834)

Description

Recombinant Human CD141 Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

THBD

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

APAEPQPGGS QCVEHDCFAL YPGPATFLNA
SQICDGLRGH LMTVRSSVAA DVISLLLNGD
GGVGRRLWI GLQLPPGCGD PKRLGPLRGF
QWVTGDNNTS YSRWARLDLN GAPLCGPLCV
AVSAAEATVP SEPIWEEQQC EVKADGFLCE
FHFPATCRPL AVEPGAAAAA VSITYGTPFA
ARGADFQALP VGSSAAVAPL GLQLMCTAPP
GAVQGHWARE APGAWDSCSVE NGGCEHACNA
IPGAPRCQCP AGAALQADGR SCTASATQSC
NDLCEHFCVP NPDQPGSYSC MCETGYRLAA
DQHRCEDVDD CILEPSPCPQ RCVNTQGGFE
CHCYPNYDLV DGECVEPVPD CFRANCEYQC
QPLNQTSYLC VCAEGFAIP HEPHRCQMFC
NQTACPADCD PNTQASCECP EGYILDDGFI
CTDIDECENG GFCSGVCHNL PGTFCICGP
DSALARHIGT DCDSGKVDGG DSGSGEPPPS
PTPGSTLTPP AVGLVSHHH HHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P07204

Entrez

7056