

Product Datasheet

Human BCAM protein (orb706857)

Description

Recombinant Human CD239 Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

BCAM

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

EVRLSVPLV EVMRGKSVIL DCTPTGTHDH
YMLEWFLTDR SGARPLASA EMQGSELQVT
MHDTRGRSPP YQLDSQRLV LAEAQVGDER
DYVCVVRAGA AGTAEATARL NVFAKPEATE
VSPNKGTLV MEDSAQEIAT CNSRNGNPAP
KITWYRNGQR LEVPVEMNPE GYMSTRVRE
ASGLLSLTST LYLRRLKDDR DASFHCAAHY
SLPEGRHGRL DSPTFHLTLH YPTEHVQFWV
GSPSTPAGWV REGDTVQLLC RGDGSPSPEY
TLFRLQDEQE EVLNVNLEGN LTLEGVTRGQ
SGTYGCRVED YDAADDVQLS KTELELVAYL
DPLELSEGKV LSLPLNSSAV VNCSVHGLPT
PALRWTKDST PLGDGPMLSL SSITFDSNGT
YVCEASLPTV PVLSTQNFT LLVQGSPELK
TAEIEPKADG SWREGDEVTL ICSARGHPDP
KLSWSQLGGS PAEIPGRQG WVSSSLTLKV
TSALSRDGI CEASNPHGNK RHVFHFGTVS
POTSOAHHHH HH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558