

Product Datasheet

Human TNFSF13 protein (orb706860)

Description

Recombinant Human CD256 Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFSF13

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

DYKDDDDKHH HHHHGPGQVQ LQKKQHSVLH
LVPINATSKD DSDVTEVMWQ PALRRGRGLQ
AQGYGVRIQD AGVYLLYSQV LFQDVTFTMG
QVVSREGQGR QETLFCIRS MPSHPDRAYN
SCYSAGVFHL HQGDILSVII PRARAKLNLS
PHGTFLGFVK LGGGGSKKQH SVLHLVPINA
TSKDDSDVTE VMWQPALRRG RGLQAQGYGV
RIQDAGVYLL YSQVLFQDVT FTMGQVVSRE
GQGRQETLFR CIRSMPSPHD RAYNSCYSAG
VFHLHQGDIL SVIIPRarak LNLSPHGTFL
GFVKLGGGGS KKQHSVLHLV PINATSKDDS
DVTEVMWQPA LRRGRGLQAQ GYGVRIQDAG
VYLLYSQVLF QDVTFTMGQV VSREGQGRQE
TLFCIRSMP SHPDRAYNSC YSAGVFHLHQ
GDILSVIIPR ARAKLNLSPH GTFLGFVKL

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date

6 months from date of receipt.