

Product Datasheet

Human DPP4 protein (orb706863)

Description

Recombinant Human CD26 Protein is

and the target gene encodi...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

DPP4

Preservatives

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

GSGSDKTHTC PPCAPELLG GPSVFLFPPK
 PKDTLMISRT PEVTCVVVDV SHEDPEVKFN
 WYVDGVEVHN AKTKPREEQY NSTYRVVSVL
 TVLHQDWLNG KEYKCKVSNK ALPAPIEKTI
 SKAKGQPREP QVYTLPPSRE EMTKNQVSLT
 CLVKGFYPSD IAVEWESNGQ PENNYKTTTP
 VLDSGDSFFL YSKLTVDKSR WQQGNVFSCS
 VMHEALHNHY TQKSLSLSPG KIEGRNKGTD
 DATADSRKTY TLTDYLKNTY RLKLYSLRWI
 SDHEYLYKQE NNILVFNAEY GNSSVFLENS
 TFDEFGHSIN DYSISPDGQF ILLEYNYVKQ
 WRHSYTASYD IYDLNKRQLI TEERIPNNTQ
 WVTWSPVGHK LAYVWNNDIY VKIEPNLPSY
 RITWTGKEDI IYNGITDWVY EEEVFSAYSA
 LWWSPNGTFL AYAQFNDTEV PLIEYSFYSD
 ESLQYPKTVR VPYPKAGAVN PTVKFFVNT
 DSLSSVTNAT SIQITAPASM LIGDHLYCDV
 TWATQERISL QWLRRIQNYS VMDICDYDES
 SGRWNCLVAR QHIEMSTTGW VGRFRPSEPH
 FTLDGNSFYK IISNEEGYRH ICYFQIDKKD
 CTFITKGTWE VIGIEALTSD YLYYISNEYK
 GMPGGRNLYK IQLSDYTKVT CLSCELNPER
 CQYYSVSFSK EAKYYQLRCS GPGLPLYTLH
 SSVNDKGLRV LEDNSALDKM LQNVQMPSKK
 LDFIILNETK FWYQMILPPH FDKSKKYPLL
 LDVYAGPCSO KADTVFERI NW ATYLASTENI

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

HMSHFQKQCF SLP