

Product Datasheet

Human TNFRSF10D protein (orb706864)

Description

Recombinant Human CD264 Protein is

system and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF10D

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

ATIPRQDEVP QQTVPQQQR RSLKEEECPA
GSHRSEYGA CNPCTEGVDY TIASNNLPSC
LLCTVCKSGQ TNKSSCTTTR DTVCQCEKGS
FQDKNSPEMC RTCRTGCPRG MVKVSNCTPR
SDIKCKNESA ASSTGKTPAA EETVTILGM
LASPYHIEGR IDPKSCDKTH TCPPCPAPEL
LGGPSVFLFP PKPKDTLMIS RTPEVTCVVV
DVSHEDPEVK FNWYVDGVEV HNAKTKPREE
QYNSTYRVVS VLTVLHQDWL NGKEYKCKVS
NKALPAPIEK TISKAKGQPR EPQVYTLPPS
RDELTKNQVS LTCLVKGFYP SDIAVEWESN
GQPENNYKTT PPVLDSDGSF FLYSKLTVDK
SRWQQGNVFS CSVMHEALHN HYTQKSLSL
PGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.