

Product Datasheet

Human TNFRSF8 protein (orb706875)

Description

Recombinant Human CD30 Protein is

and the target gene encodi...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

TNFRSF8

Preservatives

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

FPQDRPFEDT CHGNPSHYD KAVRRCCYRC
PMGLFPTQQC PQRPTDCRKQ CEPDYILDEA
DRCTACVTCS RDDLVKTPC AWNSSRVCEC
RPGMFCSTSA VNSCARCFH SVCPAGMIVK
FPGTAQKNTV CEPASPGVSP ACASPENCKE
PSSGTIPQAK PTPVSPATSS ASTMPVRGGT
RLAQEAASKL TRAPDSPSSV GRPSSDPGLS
PTQPCPEGSG DCRKQCEPDY YLDEAGRCTA
CVSCSRDDLK EKTCAWNSS RTCECRPGMI
CATSATNSCA RCVYPICAA ETVTKPQDMA
EKDTTFEAPP LGTQPCDNPT PENGEAPAST
SPTQSLLVDS QASKTLPIPT SAPVALSSTG
KIEGRDMDPK SCDKTHTCPP CPAPELLGGP
SVFLFPPKPK DTLMSRTPE VTCVVVDVSH
EDPEVKFNWY VDGVEVHNAK TKPREEQYNS
TYRVVSVLTV LHQDWLNGKE YKCKVSNKAL
PAPIEKTISK AKGQPREPQV YTLPPSRDEL
TKNOVSI TCI VKGEYPSDIA VFWESNGOPE

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Source

Human cells