

Product Datasheet

Human SLAMF6 protein (orb706882)

Description

Recombinant Human CD352 Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

SLAMF6

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

QSSLTPLMVN GILGESVTLP LEFPAGEKVN
FITWLFNETS LAFIVPHETK SPEIHVTNPK
QGKRLNFTQS YSLQLSNLKM EDTGSYRAQI
STKTSAKLSS YTLRILRQLR NIQVTNHSQL
FQNMTCELHL TCSVEDADDN VSFRWEALGN
TLSSQPNLTV SWDPRISSEQ DYTCAENAV
SNLSFSVSAQ KLCEDVKIQY TDTKIEGRMD
PKSCDKTHTC PPCPAPELLG GPSVFLFPPK
PKDTLMISRT PEVTCVVVDV SHEDPEVKFN
WYVDGVEVHN AKTKPREEQY NSTYRVVSVL
TVLHQDWLNG KEYKCKVSNK ALPAIEKTI
SKAKGQPREP QVYTLPPSRD ELTKNQVSLT
CLVKGFYPSD IAVEWESNGQ PENNYKTPPP
VLDSGDGSFFL YSKLTVDKSR WQQGNVFSCS
VMHEALHNHY TQKSLSLSPG K

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date

6 months from date of receipt.