

Product Datasheet

Human CD38 protein (orb706884)

Description

Recombinant Human CD38 Protein is

and the target gene encodi...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

CD38

Preservatives

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

VPRWRQQWSG PGTTKRFPET VLARCVKYTE
IHPEMRHVDC QSVWDAFKGA FISKHPCNIT
EEDYQPLMKL GTQTVPCNKI LLWSRIKDLA
HQFTQVQRDM FTLEDTLG Y LADDLTWCGE
FNTSKINYQS CPDWRKDCSN NPVSFVKTV
SRRFAEAACD VVHVMLNGSR SKIFDKNSTF
GSVEVHNLQP EKVQTLEAWV IHGGREDSRD
LCQDPTIKEL ESIISKRNIQ FSCKNIYRPD
KFLQCVKNPE DSSCTSEIHH HHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P28907

Entrez

952

Expiration Date

6 months from date of receipt.