

Product Datasheet

Human CD38 protein (orb706885)

Description

Recombinant Human CD38 Protein is

and the target gene encodi...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

CD38

Preservatives

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

GSGSHHHHHH VPRWRQQWSG PGTTKRFPE
VLARCVKYTE IHPEMRHVDC QSVWDAFKGA
FISKHPCNIT EEDYQPLMKL GTQTVPCNKI
LLWSRIKDLA HQFTQVQRDM FTLEDTLG
LADDLTWCGE FNTSKINYQS CPDWRKDCSN
NPVSVFWKTV SRRFAEAACD VVHVMLNGSR
SKIFDKNSTF GSVEVHNLQP EKVQTLAWV
IHGGREDSRD LCQDPTIKEL ESIISKRNQ
241FSCKNIY RPDKFLQCVK NPEDSSCTSE I

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P28907

Entrez

952

Expiration Date

6 months from date of receipt.