

## Product Datasheet

### Human ENTPD1 protein (orb706888)

## Description

Recombinant Human CD39 Protein is

and the target gene enco...

## Reactivity

Human

## Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

## Tested

ELISA, MS, SDS-PAGE, WB

## Applications

## Target

ENTPD1

## Preservatives

Liquid in a 0.2 μm filtered solution of 20mM Tris, 500mM NaCl, 10% Glycerol, pH 8.1.

## Form/Appearance

Liquid in a 0.2 μm filtered solution of 20mM Tris, 500mM NaCl, 10% Glycerol, pH 8.1.

## Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

## Note

For research use only

## Application notes

Please minimize freeze-thaw cycles.

## Protein Sequence

TQNKALPENV KYGIVLDAGS SHTSLYIKW  
PAEKENDTGV VHQVEECRVK GPGISKFVQK  
VNEIGIYLTDCMERAREVIP RSQHQETPVY  
LGATAGMRLR RMESEELADR VLDVVERSLR  
NYPFDQGAR IITGQEEGAY GWITINYLLG  
KFSQKTRWFS IVPYETNNQE TFGALDLGGA  
STQVTFVPQN QTIESPDNAL QFRLYGKDYN  
VYTHSFLCYG KDQALWQKLA KDIQVASNEI  
LRDPCFHPGY KKVNVSDLY KTPCTKRFEM  
TLPFQQFEIQ GIGNYQQCHQ SILELFNTSY  
CPYSQCAFNG IFLPPLQGDF GAFSAFYFVM  
KFLNLTSEKV SQEKVTEMMK KCAQPWEEI  
KTSYAGVKEK YLSEYCFSGT YILSLLQGY  
HFTADSWEHI HFIGKIQSD AGWTLGYMLN  
LTNMIPAEQP LSTPLSHSTY VHHHHHH

## Purity

Greater than 95% as determined by reducing SDS-PAGE.

## Source

Baculovirus

## Uniprot ID

**P49961**

## Entrez

**953**

## Expiration Date

6 months from date of receipt.