

Product Datasheet

Human SELL protein (orb706895)

Description

Recombinant Human CD62L Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

SELL

Preservatives

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 2mM CaCl₂, 2mM MgCl₂, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 2mM CaCl₂, 2mM MgCl₂, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

WTYHYSEKPM NWQRARRFCR DNYTDLVAIQ
NKAIEIYLEK TLPFSRSYYW IGIRKIGGIW
TWVGTNKS LT EEAENWGDGE PNNKKNKEDC
VEIYIKRNKD AGKWND DACH KLKAALCYTA
SCQPWSCSGH GECVEIINNY TCNCDVGYYG
PQCQFVIQCE PLEAPELGTM DCTHPLGNFS
FSSQCAFSCS EGTNLTGIEE TTCGPFGNWS
SPEPTCQVIQ CEPLSAPDLG IMNCSHPLAS
FSFTSACTFI CSEGTELIGK KKTICESSGI
WSNPSPICQK LDKSF SMIKE GDYNVDH HHHH HH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells