

Product Datasheet

Human CDON protein (orb706906)

Description

Recombinant Human CDO Protein is

and the target gene encodin...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

CDON

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

DLAPYFTSEP LSAVQKLGGP VVLHCSAQPV
TTRISWLHNG KTL DGNLEHI KIHQGTIL
SLNSSLLGYY QCLANNSIGA IVSGPATVSV
AVLGDFGSST KHVITAEKS AGFIGCRVPE
SNPKAEVRYK IRGKWLKHST ENYLILPSGN
LQILNVSLED KGSYKCAAYN PVTHQLKVEP
IGRKLLVSRP SSDDVHILHP THSQALAVLS
RSPVTLECVV SGVPAPQVYW LKDGQDIAPG
SNWRRLYSHL ATDSVDPADS GNYSCMAGNK
SGDVKYVTYM VNVLEHASIS KGLQDQIVSL
GATVHFTCDV HGNPAPNCTW FHNAQPIHPS
ARHLTAGNGL KISGVTVEDV GMYQCVDANG
IGFMHSTGRL EIENDGGFKP VIITAPVSAK
VADGDFVTLN CNASGLPVPV IRWYDSHGLI
TSHPSQVLRS KSRKSQLSRP EGLNLEPVYF
VLSQAGASSL HIQAVTQEHA GKYICEAANE
HGTTQAEASL MVVPFETNTK AETVTLPDAA
ONDDRSKRDG SETGLI SSEP VKVHPSAVES

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

NGGSPITAFK VEYKRMRTSN WLVAEDIPP