

Product Datasheet

Human MET protein (orb706908)

Description

Recombinant Human c-Met Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

MET

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

ECKEALAKSE MNVNMKYQLP NFTAETPIQN
VILHEHHIFL GATNYIYVLN EEDLQKVAEY
KTGPVLEHPD CFPCQDCSSK ANLSGGVWKD
NINMALVVDY YYDDQLISCG SVNRGTCQRH
VFPHNHTADI QSEVHCIFSP QIEEPSQCPD
CVVSALGAKV LSSVKDRFIN FVVGNTINSS
YFPDHPLHSI SVRRLKETKD GFMFLTDQSY
IDVLPEFRDS YPIKYVHAFF SNNFIYFLTV
QRETLDAQTF HTRIIRFCSI NSGLHXYMEM
PLECILTEKR KKRSTKKEVF NILQAAYVSK
PGAQLARQIG ASLNDLILFG VFAQSKPDSA
EPMDRSAMCA FPIKYVNDFF NKIVNKNVNR
CLQHFYGNH EHCNRTLLR NSSGCEARRD
EYRTEFTAL QRVDLFMGQF SEVLLTSIST
FIKGLTIAN LGTSEGRFMQ VVVSRSRGPST
PHVNFLLDSH PVSPEVIVEH TLNQNGYTLV
ITGKKITKIP LNGLGCRHFQ SCSQCLSAPP
EVOCGWCHDK CVRSFECL SG TWTOOICLPA

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

YEIHPTKSF I SGGSTITGVG KNLNSVSVPR