

Product Datasheet

Human EPHA1 protein (orb706928)

Description

Recombinant Human EPHA1 Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

EPHA1

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

KEVTLMDSK AQGELGWLLD PPKDGWSEQQ
QILNGTPLYM YQDCPMQGRD DTDHWLRSNW
IYRGEEASRV HVELQFTVRD CKSFPGGAGP
LGCKETFNLL YMESDQDVGI QLRRPLFQKV
TTVAADQSFT IRDLASGSVK LNVERCSLGR
LTRRGLYLAF HNPACVALV SVRVFYQRCP
ETLNGLAQFP DTLPGPAGLV EVAGTCLPHA
RASPRPSGAP RMHCSPDGEW LVPVGRCHCE
PGYEEGSGE ACVACPSGSY RMDMDTPHCL
TCPQQSTAES EGATICTCES GHYRAPGEGP
QVACTGPPSA PRNLSFSASG TQLSLRWEPP
ADTGGRQDVR YSVRCSQCQG TAQDGGPCQP
CGVGVHFSPP ARGLTTPAVH VNGLEPYANY
TFNVEAONGV SGLGSSGHAS TSVSISMGHA
ESLSGLSLRL VKKEPRQLEL TWAGSRPRSP
GANLTYELHV LNQDEERYQM VLEPRVLLTE
LQPDTTYIVR VRMLTPLGPG PFSPDHEFRT
SPPVSRGI TG GFHHHHHHH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558