

Product Datasheet

Human CFHR3 protein (orb706943)

Description

Recombinant Human FHR3 Protein is

and the target gene encodi...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

CFHR3

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

QVKPCDFPDI KHGGLFHENM RRPYFPVAVG
KYYSYCYDEH FETPSGSYWD YIHCTQNGWS
PAVPCLRKCY FPYLENGYNQ NYGRKFVQGN
STEVACHPGY GLPKAQTTVT CTEKGWSPTP
RCIRVRTCSK SDIEIENGFI SESSSIYILN
KEIQYKCKPG YATADGNSSG SITCLQNGWS
AQPICINSSE KCGPPPPISN GDTTSFLLKV
YVPQSRVEYQ CQPYEYELQGS NYVTCSNGEW
SEPPRCIHPC IITEENMNKN NIKLKGRSDR
KYYAKTGDTI EFMCKLGYN A NTSILSFQAV
CREGIVEYPR CEVDHHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

Q02985