

Product Datasheet

Human TNFRSF14 protein (orb706956)

Description

Recombinant Human HVEM Protein is

and the target gene encodi...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF14

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

LPSCKEDEYP VGSECCPKCS PGYRVKEACG
ELTGTVCEPC PPGTYIAHLN GLSKCLQCQM
CDPAMGLRAS RNCSTENAV CGCSPGHFCI
VQDGDHCAAC RAYATSSPGQ RVQKGGTESQ
DTLCQNCPPG TFSPNGTLEE CQHQTCSWL
VTKAGAGTSS SHWVVDIEG RMDEPKSCDK
THTCPPCPAP ELLGGPSVFL FPPKPKDTLM
ISRTPEVTCV VVDVSHEDPE VKFNWYVDGV
EVHNAKTKPR EEQYNSTYRV VSVLTVLHQD
WLNKEYKCK VSNKALPAPI EKTISKAKGQ
PREPQVYTLPSREEMTKNQ VSLTCLVKGF
YPSDIAVEWE SNGQPENNYK TTPPVLDSDG
SFFLYSKLTV DKSRWQQGNV FSCSVMEAL
HNHYTQKSLS LSPGK

Purity

Greater than 95% as determined by reducing SDS-PAGE.