

Product Datasheet

Human IDO2 protein (orb706960)

Description

Recombinant Human IDO2 Protein is

target gene encoding ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

IDO2

Preservatives

Liquid in a 0.2 μm filtered solution of 20mM Tris, 10%Glycerol, 1mM EDTA, 250mM NaCl, pH 8.0.

Form/Appearance

Liquid in a 0.2 μm filtered solution of 20mM Tris, 10%Glycerol, 1mM EDTA, 250mM NaCl, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

MEPHRPNVKT AVPLSLESYH ISEEYGFLLP
DSLKELPDHY RPWMEIANKL PQLIDAHQLQ
AHVDKMPLLS CQFLKGHREQ RLAHLVLSFL
TMGYVWQEGE AQPAEVLPRN LALPFVEVSR
NLGLPPILVH SDLVLTNWTN KDPDGFLEIG
NLETIISFPG GESLHGFIIV TALVEKEAVP
GIKALVQATN AILQPNQEAL LQALQRLRLS
IQDITKTLGQ MHDYVDPDIF YAGIRIFLSG
WKDNPAMPAG LMYEGVSQEP LKYSGGSAAG
STVLHAFDEF LGIRHSKESG DFLYRMRDYM
PPSHKAFIED IHSAPSLRDY ILSSGQDHLL
TAYNQCVQAL AELRSYHITM VTKYLITAAA
KAKHGKPNHL PGPPQALKDR GTGGTAVMSF
LKSVRDKTLE SILHPRGLEH HHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[Q6ZQW0](#)

Entrez

[169355](#)

Expiration Date

6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558