

Product Datasheet

Human IL2 protein (orb706979)

Description

Recombinant Human IL-2 Protein is produced

gene encoding ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

IL2

Preservatives

Lyophilized from a 0.2 μm filtered solution of 10mM Acetata-Na, 5% Trehaiose, pH 4.5.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of 10mM Acetata-Na, 5% Trehaiose, pH 4.5.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH2O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MPTSSSTKKT QLQLEHLLLD LQMILNGINN
YKNPKLTRML TFKFYMPKKA TELKHLQCLE
EELKPLEEVL NLAQSKNFHL RPRDLISNIN
VIVLELKGSE TTFMCEYADE TATIVEFLNR
WITFSQSIIS TLT

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

P60568

Entrez

3558

Expiration Date

6 months from date of receipt.