

Product Datasheet

Human IL23R protein (orb706985)

Description

Recombinant Human IL-23R Protein is

and the target gene enco...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

IL23R

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

GITNINCSGH IWVEPATIFK MGMNISIYCQ
AAIKNCQPRK LHFYKNGIKE RFQITRINKT
TARLWYKNFL EPHASMYCTA ECPKHFQETL
ICGKDISSGY PPDIPDEVTC VIYEYSGNMT
CTWNAGKLT Y IDTKYVVHV K SLETEEEQQY
LTSSYINIST DSLQGGKKYL VVWQAANALG
MEESKQLQIH LDDIVIPSAA VISRAETINA
TVPKTIYWD SQTIEKVSC EMRYKATTNQ
TWNVKEFDN FTYVQQSEFY LEPNIKYVFQ
VRCQETGKRY WQPWSSPFFH KTPETVPQVT
SKAFQHDTWN SGLTVASIST GHLTSDNRGD
VDDIEGRMDE PKSCDKTHC PPCAPELLG
GPSVFLFPPK PKDTLMISRT PEVTCVVVDV
SHEDPEVKFN WYVDGVEVHN AKTKPREEQY
NSTYRVVSVL TVLHQDWLNG KEYKCKVSNK
ALPAIEKTI SKAKGQPREP QVYTLPPSRE
EMTKNQVSLT CLVKGFYPSD IAVEWESNGQ
PENNYKTTTP VI DSDGSEFI YSKI TVDKSR

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558