

Product Datasheet

Human KLK2 protein (orb706994)

Description

Recombinant Human KLK2 Protein is

and the target gene encodi...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

KLK2

Preservatives

Liquid in a 0.2 μm filtered solution of 20mM Citrate, 150mM NaCl, pH 3.5.

Form/Appearance

Liquid in a 0.2 μm filtered solution of 20mM Citrate, 150mM NaCl, pH 3.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

PLIQSRIVGG WECEKHSQPW QVAVYSHGWA
HCGGVLVHPQ WVLTAHCLK KNSQVWLGRH
NLFEPEDTGQ RVPVSHSFPH PLYNMSLLKH
QSLRPDEDSS HDLMLLRLSE PAKITDVVKV
LGLPTQEPAL GTTCYASGWG SIEPEEFLRP
RSLQCVSLHL LSNDMCARAY SEKVTEFMLC
AGLWTGGKDT CGGDSGGPLV CNGVLQGITS
WGPEPCALPE KPAVYTKVVH YRKWIKDTIA
ANPVDHHHHH H

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P20151

Entrez

3817

Expiration Date

6 months from date of receipt.