



## Product Datasheet

### MAp2 (phospho-Ser136) antibody (orb7070)

|                            |                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | MAp2 (phospho-Ser136) antibody                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                                 |
| <b>Reactivity</b>          | Human, Rat                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                           |
| <b>Tested Applications</b> | FC, IF, IHC-P                                                                                                          |
| <b>Immunogen</b>           | KLH conjugated synthesised phosphopeptide derived from human MAP2 around the phosphorylation site of Ser136 (PP(p-S)P) |
| <b>Target</b>              | Phospho-MAP2 (Ser136)                                                                                                  |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                                  |
| <b>Form/Appearance</b>     | Liquid                                                                                                                 |
| <b>Concentration</b>       | 1mg/ml                                                                                                                 |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                        |
| <b>Note</b>                | For research use only                                                                                                  |
| <b>Isotype</b>             | IgG                                                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                                             |
| <b>MW</b>                  | 70/201 kDa                                                                                                             |
| <b>Uniprot ID</b>          | <b>P11137</b>                                                                                                          |
| <b>Dilution Range</b>      | IHC-P=1:100-500, IF=1:100-500, Flow-Cyt=0.2µg /Test, IHC-F=1:100-500, ELISA=1:5000-10000                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                        |



Flow cytometry  
of RSC96  
Cells(blue).Prim...



Flow cytometry  
of SH-  
SY5Y.Primary  
Antibo...