

Product Datasheet

Human SERPINE2 protein (orb707049)

Description

Recombinant Human Serpin E2 Protein is

and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg)
as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

SERPINE2

Preservatives

Liquid in a 0.2 μm filtered solution of 20mM
Tris, 150mM NaCl, pH 8.0.

Form/Appearance

Liquid in a 0.2 μm filtered solution of 20mM
Tris, 150mM NaCl, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is
recommended that the protein be aliquoted
for optimal storage. Avoid repeated freeze-
thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

SHFNPLSLEE LGSNTGIQVF NQIVKSRPHD
NIVISPHGIA SVLGMLQLGA DGRTHKQLAM
VMRYGVNGVG KILKKINKAI VSKKNKDIVT
VANAVFVKNA SEIEVPFVTR NKDVVFQCEVR
NVNFEDPASA CDSINAWVKN ETRDMIDNLL
SPDLIDGVLT RLVLVNAVYF KGLWKSFRQP
ENTKKRTFVA ADGKSYQVPM LAQLSVFRCG
STSAPNDLWY NFIELPYHGE SISMLIALPT
ESSTPLSAII PHISTKTIDS WMSIMVPKRV
QVILPKFTAV AQTDLKEPLK VLGITDMFDS
SKANFAKITR SENLHVSHIL QKAKIEVSED
GTKASAATTA ILIARSSPPW FIVDRPFLFF
IRHNPTGAVL FMGQINKPHH HHHH

Purity

Greater than 95% as determined by
reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P07093

Entrez

5270

Expiration Date

6 months from date of receipt.