

Product Datasheet

Human VSIG8 protein (orb707075)

Description

Recombinant Human VSIG8 Protein is

and the target gene encod...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

VSIG8

Preservatives

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

VRINGDGQEV LYLAEGDNVR LGCPYVLDPE
DYGPNGLDIE WMQVNSDPAH HRENVFLSYQ
DKRINHGSLP HLQQRVRFAA SDPSQYDASI
NLMNLQVSDT ATYECRVKKT TMATRKVIVT
VQARPAVPMC WTEGHMTYGN DVVLKCYASG
GSQPLSYKWA KISGHHYPYR AGSYTSQHSY
HSELSYQESF HSSINQGLNN GDLVLKDISR
ADDGLYQCTV ANNVGYSVCV VEVKVSDSRR
IGVDDIEGRM DEPKSCDKTH TCPPCPAPEL
LGGPSVFLFP PKPKDTLMIS RTPEVTCVVV
DVSHEDPEVK FNWYVDGVEV HNAKTKPREE
QYNSTYRVVS VLTVLHQDWL NGKEYKCKVS
NKALPAPIEK TISKAKGQPR EPQVYTLPPS
REEMTKNQVS LTCLVKGFYP SDIAVEWESN
GQPENNYKTT PPVLDSGGSF FLYSKLTVDK
SRWQQGNVFS CSVMHEALHN HYTQKSLSLS
PGK

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558