

Product Datasheet

Mouse SERPINF2 protein (orb707083)

Description

Recombinant Mouse Alpha-2-antiplasmin

expression system and the tar...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

SERPINF2

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, 20% glycerol, 1mM DTT, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, 20% glycerol, 1mM DTT, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

VDLPGQPVS EQAQQKLPLP ALFKLDNQDF
GDHATLKRSP GHCKSVPTAE ETRRLAQAMM
AFTTDLFSLV AQTSTSSNLV LSPLSVALAL
SHLALGAQNQ TLHSLHRVLH MNTGSCPLHL
LSHFYQNLGP GTIRLAARIY LQKGFPIKDD
FLEQSERLFG AKPVKLTGKQ EEDLANINQW
VKEATEGKIE DFLSELPDST VLLLLNAIHF
HGFWRTKFDP SLTQKDFHHL DERFTVSVDM
MHA VSYPLRW FLLEQPEIQV AHFPFKNNMS
FVVVMPTYFE WNVSEVLANL TWDTLYHPSL
QERPTKVWLP KLHLQQQLDL VATLSQLGLQ
ELFQGPDLRG ISEQNLVVSS VQHQSTMELS
EAGVEAAAAT SVAMNRM SLS SFTVNRPF LF
FIMEDTIGVP LFGVSVRNPN PSALPQLQEQ
RDSPDNRLIG QNDKADF HGG KTFGPDLKLA
PRMEEDYPQF SSPKVDH HHH HH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q61247](#)

Entrez

[18816](#)

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558