

Product Datasheet

Mouse BTNL2 protein (orb707089)

Description

Recombinant Mouse BTNL2 Protein is

and the target gene encod...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

BTNL2

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

DDFRVVGPNL PILAKVGEDA LLTCQLLPKR
TTAHMEVRWY RSDPDMPVIM YRDGA EVTGL
PMEGYGGRAE WMEDSTEEGS VALKIRQVQP
SDDGQYWCRF QEGDYWRETS VLLQVAALGS
SPNIHVEGLG EGEVQLVCTS RGWFPEPEVH
WEGIWGEKLM SFSENHPGE DGLFYVEDTL
MVRNDSVETI SCFIYSHGLR ETQEATIALS
ERLQTELASV SVIGHSQPSP VQVGENIELT
CHLSPQTDAQ NLEVRWLRSR YYPVHVYAN
GTHVAGEQMV EYKGR TSLVT DAIHEGKLTL
QIHNARTSDE GQYRCLFGKD GYVQEARVDV
QVMAVGSTPR ITREVLKDGG MQLRCTSDGW
FPRPHVQWRD RDGKTMP SFS EAFQQGSQEL
FQVETLLLVT NGSMVNV TCS ISLPLGQ EKT
ARFPLSHHHH HH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date

6 months from date of receipt.