

Product Datasheet

Mouse THBD protein (orb707106)

Description

Recombinant Mouse CD141 Protein is

and the target gene encod...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

THBD

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0 .

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0 .

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

LSALAKLQPT GSQCVEHECF ALFQGPATFL
DASQACQRLQ GHLMTVRSSV AADVISLLLS
QSSMDLGPWI GLQLPQGCDD PVHLGPLRGF
QWVTGDNHTS YSRWARPNDQ TAPLCGPLCV
TVSTATEAAP GEPAWEEKPC ETETQGFLCE
FYFTASCRPL TVNTRDPEAA HISSTYNTPF
GVSGADFQTL PVGSSAAVEP LGLELVCRAP
PGTSEGHAW EATGAWNCSV ENGGCEYLCN
RSTNEPRCLC PRDMDLQADG RSCARPVVQS
CNELCEHFCV SNAEVPGSYS CMCETGYQLA
ADGHRCEDVD DCKQGPNPCP QLCVNTKGGF
ECFCYDGYEL VDGECEVELLD PCFGSNCEFQ
CQPVSPDYR CICAPGFAPK PDEPHKCEMF
CNETSCPADC DPNSPTVCEC PEGFILDEGS
VCTDIDCSQ GECFTSECRN FPGSYECICG
PDTALAGQIS KDCDPIPVRE DTKEEEGSGE
PPVSPTPGSP TGPPSARPVH SHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P15306

Entrez

21824