

Product Datasheet

Mouse ENTPD2 protein (orb707153)

Description

Recombinant Mouse CD39L1 Protein is

and the target gene enco...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

ENTPD2

Preservatives

Liquid in a 0.2 μm filtered solution of 50mMTris, 10mMNaCl, 150mMNaCl, pH 7.5, 10%glycerol .

Form/Appearance

Liquid in a 0.2 μm filtered solution of 50mMTris, 10mMNaCl, 150mMNaCl, pH 7.5, 10%glycerol .

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

CVPTQDVREP PALKYGIVLD AGSSHTSMFV
YKWPADKEND TGIVGQHSSC DVRGGGISSY
ANDPSRAGQS LVECLEQALR DVPKDRYAST
PLYLGATAGM RLLNLTSP EA TAKVLEAVTQ
TLTRYPFDFR GARILSGQDE GVFGWVTANY
LLENFIKYGW VGRWIRPRKG TLGAMD LGGA
STQITFETTS PSEDPDNEVH LRLYGQHYRV
YTHSFLCYGR DQVLQRLLAS ALQIHRFHPC
WPKGYSTQVL LREYQSPCT MGQRPQTFNS
SATVSLSGTS NAALCRDLVS GLFNISSCPF
SQCSFNGVFQ PPVAGNFIAF SAFYYTVDFL
KTMGLPVG T LKQLEDATET TCNQTWAE LQ
ARVPGQQTRL PDYCAVAMFI HQLLSRGYRF
DERSFRGVVF EKKAADTAVG WALGYMLNLT
NLIPADLPGL RKGTHFSHHH HHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[O55026](#)

Entrez

[12496](#)