

Product Datasheet

Mouse CD62P protein (orb707162)

Description

Recombinant Mouse CD62P Protein is

and the target gene encod...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

SELP

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

WTYNYSTKAY SWNNSRVFCR RHFTDLVAIQ
NKNEIAHLND VIPFFNSYYW IGIRKINNKW
TWVGTNKTLT EEAENWADNE PNNKKNQDC
VEIYKNSA PGKWNDPCF KRKRALCYTA
SCQDMSCSNQ GECIETIGSY TCSCYPGFY
PECEYVKECG KVNIPQHVLN NCSHPLGEFS
FNSQCTFSCA EGYELDGPGE LQCLASGIWT
NNPPKCDVAV QQSLEAPPHG TMACMHPIAA
FAYDSSCKFE CQPGYRARGS NTLHCTGSGQ
WSEPLPTCEA IACEPPEIPI HGSMDCVPST
GTFGYNSSCT FLCAEGFVLK GNDAIQCADS
GQWTAPAPFC EALQCPEFPV PSKAQVNCSD
PFGTLTYQSV CSFSCDEGSL LVGASVIRCL
ATGHWNGAPP ECQAVSCAPM LSPENGSMTC
VQPLGNSTYK STCQFMCDEG FYLSGPERLD
CSPSGHWTGT PPTCEAIKCP GIFAPEQGNL
DCSHVHGEFG VGSICHFSCN EDFELLGSEN
VECTVSGRWS APPPTCKGIT SI PAPAVERCP

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Purity

Greater than 95% as determined by reducing