

## Product Datasheet

### Mouse CHIL3 protein (orb707177)

## Description

Recombinant Mouse CHIL3 Protein is

and the target gene encod...

## Reactivity

Mouse

## Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

## Tested Applications

ELISA, MS, SDS-PAGE, WB

## Target

CHIL3

## Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

## Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

## Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

## Note

For research use only

## Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## Protein Sequence

YQLMCYYTSW AKDRPIEGSF KPGNIDPCLC  
THLIYAFAGM QNNEITYTHE QDLRDYEALN  
GLKDKNTEK TLLAIGGWKF GPASFAMVS  
TPQNRQIFIQ SVIRFLRQYN FDGLNLDWQY  
PGSRGSPPKD KHLFSLVKE MRKAFEEESV  
EKDIPRLLLT STGAGIIDVI KSGYKPELS  
QSLDYIQVMT YDLHDPKDG YTGENSEPLYKS  
PYDIGKSADL NVDSIISYWK DHGAASEKLI  
VGFPAYGHTF ILSDPSKTGI GAPTISTGPP  
GKYTDESGLL AYYEVCTFLN EGATEVWDAP  
QEVPIYQGN EWVGYDNVRS FKLKAQWLKD  
NNLGGAVVWP LDMDDFSGSF CHQRHFPLTS  
TLKGDNIHS ASCKGPYVDH HHHHH

## Purity

Greater than 95% as determined by reducing SDS-PAGE.

## Source

Human cells