

Product Datasheet

Mouse PPIA protein (orb707185)

Description

Recombinant Mouse Cyclophilin A Protein is

target gene ...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

PPIA

Preservatives

Liquid in a 0.2 μm filtered solution of PBS, 10%glycerol, pH 7.4.

Form/Appearance

Liquid in a 0.2 μm filtered solution of PBS, 10%glycerol, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

MVNPTVFFDI TADDEPLGRV SFELFADKVP
KTAENFRALS TGEKGFGYKG SSFHRIIPGF
MCQGGDFTRH NGTGGRSIYG EKFEDEFIL
KHTGPGILSM ANAGPNTNGS QFFICTAKTE
WLDGKHVVFG KVKEGMNIVE AMERFGSRNG
KTSKITISD CGQL

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[P17742](#)

Entrez

[268373](#)

Expiration Date

6 months from date of receipt.