

Product Datasheet

Mouse FSTL1 protein (orb707204)

Description

Recombinant Mouse FSTL1 Protein is

and the target gene encod...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

FSTL1

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

EEPRSKSKI CANVFCGAGR ECAVTEKGEP
TCLCIEQCKP HKRPVCGSNG KTYLNHCELH
RDACLTGSKI QVDYDGHCKE KKSASPSASP
VVCYQANRDE LRRRLIQWLE AEIIPDGWFS
KGSNYSEILD KYFKSFDNGD SHLDSSEFLK
FVEQNETAIN ITTYADQENN KLLRSLCVDA
LIELSDENAD WKLSFQEFLK CLNPSFNPPE
KKCALEDETY ADGAETEVCDC NRCVCSCGHW
VCTAMTCDGK NQKGVQTHTE EEKTGYVQEL
QKHQGTAEKT KKVNTKEIVD DIEGRMDEPK
SCDKTHTCPP CPAPELLGGP SVFLFPPKPK
DTLMISRTPE VTCVVVDVSH EDPEVKFNWY
VDGVEVHNAK TKPREEQYNS TYRVVSVLTV
LHQDWLNGKE YKCKVSNKAL PAPIEKTISK
AKGQPREPQV YTLPPSREEM TKNQVSLTCL
VKGFPYSDIA VEWESNGQPE NNYKTTTPVL
DSDGSFFLYS KLTVDKSRWQ QGNVFSCSVM
HEALHNHYTO KSI SI SPGK

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558