

Product Datasheet

Mouse HPX protein (orb707209)

Description

Recombinant Mouse Hemopexin Protein is

and the target gene e...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

HPX

Preservatives

Liquid in a 0.2 μm filtered solution of 20mM MES, 150mM NaCl, pH 5.5.

Form/Appearance

Liquid in a 0.2 μm filtered solution of 20mM MES, 150mM NaCl, pH 5.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

SPLPTANGRV AEVENGTKPD SDVPEHCLDT
WSFDAATMDH NGTMLFFKGE FVWRGHSGTR
ELISARWKNP ITSVDAAFRG PDSVFLIKED
KVWVYPPEKK ENGYPKLFQE EFPGIPYPPD
AAVECHRGEC QSEGVLFQFQ NRKWFWDFA
RTQKERSWST VGNCTAALRW LERYYCFQGN
KFLRFNPVTG EVPPRYPLDA RDYFVSCPGR
GHGRPRNGTA HGNSTHPMHS RCSPDPGLTA
LLSDHRGATY AFTGSHYWRL DSSRDGWHWS
PIAHHWPQGP STVDAAFSWD DKVYLIQGTQ
VYVFLTKGGN NLVSGYPKRL EKELGSPPGI
SLETIDAAFS CPGSSRLYVS SGRRLWWLDL
KSGAQATWTE VSWPHEKVDG ALCLDKSLGP
NTCSSNGSSL YFIHGPNLYC YSSIDKLNA
KSLPQPQKVN SILGCSQHHH HHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q91X72](#)

Entrez

[15458](#)

Expiration Date

6 months from date of receipt.