

Product Datasheet

Mouse MMP9 protein (orb707254)

Description

Recombinant Mouse MMP9 Protein is

and the target gene encodi...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

MMP9

Preservatives

Liquid in a 0.2 μm filtered solution of 20mM Tris, 150mM NaCl, pH 7.5.

Form/Appearance

Liquid in a 0.2 μm filtered solution of 20mM Tris, 150mM NaCl, pH 7.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

APYQRQPTFV VFPKDLKTSN LTDTQLAEAY
LYRYGYTRAA QMMGEKQSLR PALLMLQKQL
SLPQTGELDS QTLKAIRTPR CGVPDVGRFQ
TFKGLKWDHH NITYWIQNYE EDLPRDMIDD
AFARAFVWVG EVAPLTFTRV YGPEADIVIQ
FGVAEHGDGY PFDGKDGLLA HAFPPGAGVQ
GDAHFDDDEL WSLGKGVVIP TTYGNSNGAP
CHFPFTFEGR SYSACTTDGR NDGTPWCSTT
ADYDKDGKFG FCPSELYTE HNGEGKPCV
FPFIFEGRSY SACTTKGRSD GYRWCATTAN
YDQDKLYGFC PTRVDATVVG GNSAGELCVF
PFVFLGKQYS SCTSDGRRDG RLWCATTSNF
DTDKKWGFCP DQGYSFLVA AHEFGHALGL
DHSSVPEALM YPLYSYLEGF PLNKDDIDGI
QYLYGRGSKP DPRPPATTTT EPQPTAPPTM
CPTIPPTAYP TVGPTVGPTG APSPGPTSSP
SPGPTGAPSP GPTAPPTAGS SEASTESLSP
ADNPCNVDVF DAIAEIQGAL HFFKDGWYWK
FLNHRGSPLQ GPFLTARTWP ALPATLDSAF
EDPQTKRVFF FSGRQMWWYT GKTVLGPRSL
DKLGLGPEVT HVSGLLPRL GKALLFSKGR
VWRFDLKSQK VDPQSVIRVD KEFSGVPWNS
HDIFQYQDKA YFCHGKFFWR VSFQNEVNKV
DHEVNQVDDV GYVTDLLQC PGGGGSHHHH
HHHHHH

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558