

Product Datasheet

Mouse PLA2G1B protein (orb707266)

Description

Recombinant Mouse PLA2G1B Protein is

and the target gene enc...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

PLA2G1B

Preservatives

Liquid in a 0.2 μm filtered solution of 20mM HEPES, 150mM NaCl, pH 7.0.

Form/Appearance

Liquid in a 0.2 μm filtered solution of 20mM HEPES, 150mM NaCl, pH 7.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

AHSISPRAVW QFRNMIKCTI PGSDPLKDYN
NYGCYCGLGG WGTPVDDLDR CCQTHDHCYS
QAKKLESCKF LIDNPYTNTY SYSCSGSEIT
CSAKNNKCED FICNCDREAA ICFSKVPYNK
EYKNLDTGKF CHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

[Q9Z0Y2](#)

Entrez

[18778](#)

Expiration Date

6 months from date of receipt.