

Product Datasheet

Mouse S100A15A protein (orb707276)

Description

Recombinant Mouse S100A15A Protein is

the target gene encod...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg)
as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

S100A15A

Preservatives

Lyophilized from a 0.2 μm filtered solution
of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution
of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -
20°C, though stable at room temperature
for 3 weeks. Reconstituted protein solution
can be stored at 2-8°C for 2-7 days. Aliquots
of reconstituted samples are stable at -20°C
for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do
not mix by vortex or pipetting. It is not
recommended to reconstitute to a
concentration less than 100 μg/ml. Dissolve
the lyophilized protein in ddH₂O. Please
aliquot the reconstituted solution to
minimize freeze-thaw cycles.

Protein Sequence

MPDTPVEDSL FQIIHCFHHY AAREGDKETL
SLEELKALLL DSVPRFMDTL GRRQPYITE
LFRAADKNKD NQICFDEFLY ILGKLVKDYH
LQFHRQLCAH YCTEHSLY

Purity

Greater than 95% as determined by
reducing SDS-PAGE.

Source

E. coli

Uniprot ID
Q6S5I3
Entrez
381493
Expiration Date

6 months from date of receipt.