

Product Datasheet

Mouse TGF beta Receptor 2 protein (orb707290)

Description

Recombinant Mouse TGFBR2 Protein is

system and the target gene enco...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

TGFBR2

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

IPPHVPKSVN SDVMASDNGG
AVKLPQLCKF CDVRLSTCDN
QKSCMSNCSI TAICEKPHEV
CVAVWRKNDK NITLETVCHD
PKLTYHGFTL EDAASPKCVM KEKKRAGETF
FMCACNMEEC NDYIIFSEY TTSSPDVDDI
EGRMDEPKSC DKTHTCPPCP
APELLGGPSV FLFPPKPKDT LMISRTPEVT
CVVVDVSHED PEVKFNWYVD
GVEVHNAKTK PREEQYNSTY
RVVSVLTVLH QDWLNGKEYK
CKVSNKALPA PIEKTISKAK GQPREPQVYT
LPPSREEMTK NQVSLTCLVK GFYPSDIAVE
WESNGQPENN YKTTTPVLDS
DGSFFLYSKL TVDKSRWQQG
NVFSCSVME ALHNHYTQKS LSLSPGK

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558