

Product Datasheet

Mouse TIMP2 protein (orb707297)

Description

Recombinant Mouse TIMP2 Protein is

and the target gene encod...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

TIMP2

Preservatives

Liquid in a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Liquid in a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Please minimize freeze-thaw cycles.

Protein Sequence

CSCSPVHPQQ AFCNADVIR AKAVSEKEVD
SGNDIYGNPI KRIQYEIKQI KMFKGPDKDI
EFIYTAPSSA VCGVSLDVGG KKEYLIAGKA
EGDGKMHITL CDFIVPWDITL SITQKKSLNH
RYQMGCECKI TRCPMIPCYI SSPDECLWMD
WVTEKSINGH QAKFFACIKR SDGSCAWYRG
AAPPKQEFLL IEDPHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

P25785

Entrez

21858

Expiration Date

6 months from date of receipt.