


Product Datasheet

Talin (phospho-Ser425) antibody (orb7073)

Description	Talin (phospho-Ser425) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Talin around the phosphorylation site of Ser425 (SV(p-S)PK)
Target	Phospho-Talin (Ser425)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	280 kDa
Uniprot ID	Q9Y490
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Paraffin-embedded human endometrium carc...