

Product Datasheet

Mouse TYRO3 protein (orb707302)

Description

Recombinant Mouse TYRO3 Protein is

and the target gene encod...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications

Target

TYRO3

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

AGLKLMGAPV KMTVSQGPV KLNCSVEGME
DPDIHWMKDG TVVQNASQVS ISISEHSWIG
LLSLKSVRS DAGLYWCQVK DGEETKISQS
VWLTVEGVPP FTVEPKDLAV PPNAPFQLSC
EAVGPPEPVT IYWWRGTLKV GGPAPSPSVL
NVTGVTQRTE FSCEARNIKG LATSRAIVR
LQAPPAAPFN TTVTTISSYN ASVAWVPGAD
GLALLHSCTV QVAHAPGEWE ALAVVVPVPP
FTCLLRNLAP ATNYSRLVRC ANALGPSYGY
DWVPFQTKGL APARAPQNFH AIRTDSGLIL
EWEEVIPEDP GEGPLGPYKL SWVQENGTDQ
ELMVEGTRAN LTDWDPQKDL ILRVCSNAI
GDGPWSQPLV VSSHDHAGRQ GPPHSRTSVD
DIEGRMDEPK SCDKTHTCPP CPAPELLGGP
SVFLFPPKPK DTLMSRTPE VTCVVVDVSH
EDPEVKFNWY VDGVEVHNAK TKPREEQYNS
TYRVVSVLTV LHQDWLNGKE YKCKVSNKAL
PAPIEKTISK AKGOPREPOV YTI PPSREEM

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558