

Product Datasheet

Mouse VSIG4 protein (orb707306)

Description

Recombinant Mouse VSIG4 Protein is

and the target gene encod...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

VSIG4

Preservatives

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

HPTLKTPEVS TGTWKGDV KI QCIYDPLRGY
RQVLVKWLVR HGSDSVTIFL RDSTGDHIQQ
AKYRGRLKVS HKVPGDVSLQ INTLQMDDRN
HYTCEVTWQT PDGNQVIRDK IIELRVRKYN
PPRINTEAPT TLHSSLEATT IMSSTSDLTT
NGTGKLEETI AGSGRNLP HH HHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID

F6TUL9

Entrez

278180

Expiration Date

6 months from date of receipt.