


Product Datasheet

ATP Microplate Assay Kit (orb707328)

Description	Detection and Quantification of ATP Content.
Reactivity	All
Range	0.02 mmol/L - 2.5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 660 nm
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	56-65-5
Expiration Date	Please enquire.



Representative
Standard
Curve