


Product Datasheet

Glycogen Branching Enzyme Microplate Assay Kit (orb707332)

Description	Detection and Quantification of Glycogen Branching Enzyme Activity.
Reactivity	All
Range	0.05 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 660 nm
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-97-2
Expiration Date	Please enquire.



Representative
Standard
Curve